

M52110S-WH

OP-AIR PRO OXYGEN

RESPIRATORY FFP2 MASK with Elastic Head loops

Category III PPE and Class I, Medical Device

FEATURES & BENEFITS

- FFP2 NR D efficiency class filtering mask designed for protect from solid and liquid particles
- Ultrasonic welding

Recommendations :

- Use in a dusty environment
- Protection for those at risk of infectious diseases (According to the national recommendations in force)
- Duration of use: 8 hours maximum (continuous)




Article M52110S-WH

Color White

TECHNICAL DIESCRIPTION

Product Name: **KOLMI – Op-Air Pro Oxygen – Mask FFP2 NR D – Type IIR**
 Product Type: Single-use, non-sterile
 Int./Ext. Layers : Spunbonded polypropylene
 Filter: Melt blown polypropylene
 Link: Synthetic Elastic Head Loops
 Nasal bar: Polypropylene and metal
 Size: Large
 Unitary Weight: 6,5 g (± 10%)
 Country of Origin: France



EN 149:2001+A1:2009

Performances	Required Level	Laboratory Name	Report Number	Date of Report	Results
Paraffin oil penetration	< 6 % after 120 mg of exposure	APAVE	19.0029	22/01/2019	Pass
NaCl penetration	< 6 % after 120 mg of exposure				Pass
Facial leaks	46 results ≤ 11% 8 averages out of 10 ≤ 8%				Pass
Respiratory resistance Inhalation 30 l/min	≤ 0,7 mbars				Pass
Respiratory resistance Inhalation 95 l/min	≤ 2,4 mbars				Pass
Respiratory resistance exhalation 160 l/min	≤ 3 mbars				Pass
Carbon dioxide content	< 1,0 %				Pass
Flammability	Must not burn or continue to burn for more than 5 seconds after removing the flame				Pass

Protection (D): protection against solid and liquid aerosols, associated with superior resistance to clogging tested with dolomite dust.

* Average test results (Reception status + Simulated port processing)

Annual monitoring according to module D - Reg (EU) 2016/425 performed by APAVE

NF EN 14683:2019

TEST	STANDARD	REQUIRED LEVEL	LABORATORY	REPORT NUMBER	DATE OF REPORT	RESULTS
BACTERIAL FILTRATION EFFICIENCY (BFE %)	EN 14683:2019	≥ 98% (Type IIR)	CENTEXBEL	20.00718.02	20/03/2020	99.9% min
DELTA P	EN 14683:2019	< 60 Pa/cm ²	CENTEXBEL	20.00718.02	20/03/2020	46,3 Pa/cm ² max
SPLASH	ISO 22609:2004	≥ 16 kPa	CENTEXBEL	20.00718.02	20/03/2020	PASS
Cytotoxicity	ISO 10993-5	Absence de cytotoxicité	NELSON	1030830-S01	20/03/2018	Absence of cytotoxicity
Irritation cutaneous test	ISO 10993-10	Non irritant	NAMSA	231037	21/06/2017	Non irritant
Sensibilization	ISO 10993-10	Non sensibilisant	NAMSA	231038	19/06/2017	No sensibilization
Microbial sanitation	NF 11737:2018	≤30 cfu/g	MICROSEPT	1812011-15	14/04/2020	PASS
Duration of use BFE + Delta P	EN 14683:2014	≥ 98% (Type IIR) < 49 Pa/cm ²	NELSON	961620-S01	28/04/2017	8h : > 99.9% 8h : < 37,00 Pa/cm ²

PRECAUTIONS FOR USE

The device should only be used on healthy skin.
Their re-use or prolonged use can produce infection or cross-contamination.
After use, comply with the national regulations in force for the disposal of the device.

CERTIFICATION & STANDARDS

Meets the requirements of Regulation (EU) 2017/745 and Regulation (EU) 2016/425.
Complies with the applicable harmonized standards EN 14683 and EN 149.
Manufactured under an ISO 13485 and ISO 9001 certified system.

STORAGE

Normal conditions of conservation and storage: must not be exposed to humidity and sun, must be stored at a temperature between 5°C and 40°C.
Product life: 5 years.

LOGISTICS



Carton Specifications

Dispenser Box Specifications

Article	Size mm	Unit weight (kg)	QTY/ Palette	Size mm	Cdts	QTY
M52110S-WH	355 x 235 x 330	1,9	45 cartons (4 layers of 10 cartons + 1 layer of 5 cartons)	310 x 110 x 85	Individual Packet	4 dispenser boxes of 50 units

