



# SafeMask®

●●●●● FFP2 Masks



## Protection that lets you breathe easy



# Quality ... made in Europe

3 SIZES  
FOR OPTIMAL  
PROTECTION

## Individual fit

Medicom FFP2 masks are available in three sizes. A special size guide helps each wearer select the correct size for optimal protection.



## Exceptionally breathable


High-quality raw materials offer exceptional breathability. Non-woven spunbond materials from leading manufacturers ensure that Medicom FFP2 masks provide reliable protection and maximum comfort, even during extended wear.

## Superior comfort


Ultra-soft, latex-free inner layer makes this mask ideal for those with sensitive skin.

# One mask, twice the protection

Medicom FFP2 masks are unique: their high bacterial filtration efficiency (BFE 99.9%) means they protect the wearer against solid and liquid aerosols, while their Type IIR classification means they also protect wearers against contamination from splashes and sprays.



Two ultrasonically sealed elastic headbands create an optimal fit and a reliable seal around the face.




Spunbond and meltblown layers ensure maximum protection and unrivaled comfort.

Max. wear time: 8 hours.

## Sizing for all

To maximize efficacy, Medicom FFP2 masks are available in three sizes.



# How to order

+31 33 286 6460  
@ cs@medicom-eu.com



## SafeMask® Respiratory Face Mask EU Standard FFP2/Type IIR



### Medicom FFP2 Mask

PPE/MEDICAL FFP2 NR D  
EN 149, Type IIR EN 14683

Colour: white

Content: 50 pcs.

Size: REF:

(S) Small 2292-WH

(M) Medium 2092-WH

(L) Large 2192-WH

Content: 50 pcs. (single packed)

Size: REF:

(S) Small 2292S-WH

(M) Medium 2092S-WH

(L) Large 2192S-WH