

## 1129 SAFETOUCH® EVERSTRONG™ VINYL POWDER-FREE Class I, Medical Device and Category II, Personal Protective Equipment

### FEATURES & BENEFITS

- Suitable elasticity
- Good dexterity
- Its smooth appearance gives it correct touch and grip of objects
- Watertight and airtight

### Recommendation:


- Glove for professional use
- In the healthcare and industrial environment, this glove is used to make barrier against microorganisms between the wearer of the glove and the examined people or objects / materials handled.
- Powder-free, it reduces cross contamination and limits the risks of irritation.
- For food handling / preparations (HACCP)



### PRODUCT DESCRIPTION

Product name:	<b>SafeTouch® EverStrong™</b>
Product type:	Single use, non-sterile
Material :	Vinyl
Colour :	Transparent
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Smooth
Powdered :	No
Length :	240 mm
Palm thickness :	0.08 ±0,02 mm
Finger thickness :	0.09 ±0,02 mm
Unit weight (M):	5,3±0,2 g
Country of origin :	China
AQL :	<1,5
Strength at break :	
✓ Before aging	≥ 3,6 N
✓ After aging	≥ 3,6 N

Item	Size	
1129-A	T 6	S
1129-B	T 7	M
1129-C	T 8	L
1129-D	T 9	XL



### INDUSTRIES

Food / Non residential catering / Community / Administration / Health / Medical / Dentist / Tattoo / Sanitary / Maintenance / Cleaning / Lab / Pharmaceutical / Research

## CERTIFICATION

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC  
Comply with the EN 455-1/2/3/4 Medical gloves for single use

Manufactured under ISO 13485 certified quality system

In accordance with requirements of personal protective equipment Dir.89/686/EEC

Comply with the EN 420 and EN 374-2:2003 : Protection against microorganisms.

In accordance with requirements of Reg.(CE) 1935/2004 and Reg.(EU) 10/2011. This glove is suitable for food contact with all food\*.

\* Restriction : It is not recommended to put vinyl powdered gloves in contact more than 30 minutes with foods containing fats or oils (see Annex 1).

## STORAGE CONDITIONS

Normal conditions of conversation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

Shelf life of the product: 5 years

## LOGISTIC INFORMATION



### Specifications outer case




### Specifications inner case

### Colour code inner and outer

Item	Size mm	Gross weight kg	QTY/Pallet	Size mm	QTY
1129-A	320 x 250 x 240	5,18	60	230 x 120 x 70	10 x 100
1129-B	320 x 250 x 240	6,03	60	230 x 120 x 70	10 x 100
1129-C	320 x 250 x 240	6,36	60	230 x 120 x 70	10 x 100
1129-D	320 x 250 x 240	6,90	60	230 x 120 x 70	10 x 100

Colour







### EAN Codes

Item	Inner case	Outer case
1129-A	3 700 381 501 963	3 700 381 500 515
1129-B	3 700 381 501 970	3 700 381 500 522
1129-C	3 700 381 501 987	3 700 381 500 539
1129-D	3 700 381 501 994	3 700 381 500 546