

1163

LATEX GLOVE SAFETOUCH® CONNECT™ REJUVENATE POWDER-FREE Enriched in Lanolin and Vitamin E Class I, Medical Device

FEATURES & BENEFITS

- Great elasticity
- Good dexterity
- Latex is the material which allows the best contact sensibility
- Its micro-textured appearance gives it correct touch and grip of objects
- Enriched in Lanolin and Vitamin E, this glove gives the feeling of a better hydrated skin, softer and smoother.
- Contraindication: Latex is a natural material that contains protein which might cause allergies



Recommendation:

- In the healthcare environment, this glove is used to make barrier against microorganisms between the wearer of the glove and the examined people or objects / materials handled.
- Powder-free, it reduces airborne cross contamination and limits the risk of irritation.

PRODUCT DESCRIPTION

Product name :	SafeTouch® ADVANCED™ Rejuvenate
Product type :	Single use, non-sterile
Material :	Latex
Colour :	Natural
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Micro-textured
Powdered :	No
Length :	240 mm
Palm thickness :	0,14±0,03 mm
Finger thickness :	0,15±0,03 mm
Integrated lotion:	Lanolin and Vitamin E
Unit weight (M):	6 g
Country of origin :	Malaysia
AQL :	< 1,5
Tensile Strength :	
✓ Before aging	≥21 MPa
✓ After aging	≥16 MPa
Elongation:	
✓ Before aging	≥500%
✓ After aging	≥500%

Item	Size	
1163-A	T 6	XS
1163-B	T 7	S
1163-C	T 8	M
1163-D	T 9	L



INDUSTRIES

Health / Medical / Lab / Pharmaceutical / Research.

CERTIFICATION

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC.

Comply with the EN 455-1/2/3/4 Medical gloves for single use.

Manufactured under ISO 13485 certified quality system.

Additional cosmetic substances (Lanolin & Vitamin E) are in accordance with the Reg.(EC) 1223/2009 and its amendments concerning cosmetics substances.

STORAGE CONDITIONS

Normal conditions of conversation & storage: not to be exposed to moisture and sunlight, must be stored at a temperature between 5°C and 40°C

Shelf life of the product : 5 years

LOGISTIC INFORMATION







Specifications outer case



Specifications inner case

Colour code inner and outer

Item	Size mm	Gross weight kg	QTY/ Pallet	Size mm	QTY	Colour
1163-A	370 x 255 x 250	6,3	63	240 x 120 x 70	10 x 100	
1163-B	370 x 255 x 250	6,8	63	240 x 120 x 70	10 x 100	
1163-C	370 x 255 x 250	7,1	63	240 x 120 x 70	10 x 100	
1163-D	370 x 255 x 250	8,3	63	240 x 120 x 70	10 x 100	

EAN Codes

Item	Inner case	Outer case
1163-A	8 718 007 902 167	8 718 007 002 164
1163-B	8 718 007 902 174	8 718 007 002 171
1163-C	8 718 007 902 181	8 718 007 002 188
1163-D	8 718 007 902 198	8 718 007 002 195