

1172 SAFETOUCH® ADVANCED™ EXTEND NITRILE POWDER-FREE Class I, Medical Device and Category II, Personal Protective Equipment

FEATURES & BENEFITS

- A touch of elegance with this pink color
- Its elasticity and flexibility provide excellent dexterity.
- Ergonomically, it's as comfortable as a latex glove
- Its textured appearance improve grip of objects and facilitates their manipulation.
- Watertight and airtight

Recommendation:

- Glove for professional use
- In the healthcare and industrial environment, this glove is used to make barrier against microorganisms between the wearer of the glove and the examined people or objects / materials handled.
- Powder-free, it reduces cross contamination and limits the risk of irritation.



PRODUCT DESCRIPTION

Product name :	SafeTouch® Advanced™ Extend
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	Pink
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Textured
Powdered :	No
Length :	240 mm
Palm thickness :	0,09 mm
Finger thickness :	0,11 mm
Unit weight (M):	4,5 ± 0,2 g
Country of origin :	Malaysia
AQL :	< 1,5
Tensile Strength :	
✓ Before aging	≥14MPa
✓ After aging	≥14MPa
Elongation:	
✓ Before aging	≥ 500%
✓ After aging	≥ 400%

Item	Size	
1172-A	T 6	XS
1172-B	T 7	S
1172-C	T 8	M
1172-D	T 9	L



INDUSTRIES

Health / Medical / Dentist / Tattoo / Community / Administration / Sanitary / Hygiene / Maintenance / Cleaning / Lab / Pharmaceutical / Research

CERTIFICATION

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC
 Comply with the EN 455-1/2/3/4 Medical gloves for single use
 Manufactured under ISO 13485 and ISO 9001 certified quality system
 In accordance with requirements of personal protective equipment Dir.89/686/EEC
 Comply with the EN 420 and EN 374-2:2003 : Protection against microorganisms.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored at a temperature between 5°C and 40°C
 Shelf life of the product: 5 years

LOGISTIC INFORMATION







Specifications outer case



Specifications inner case

Colour code inner and outer

Item	Size mm	Gross weight kg	QTY/pallet	Size mm	QTY	Colour
1172-A	365 x 260 x 250	4,8	63	240 x 120 x 70	10 x 100	
1172-B	365 x 260 x 250	5,2	63	240 x 120 x 70	10 x 100	
1172-C	365 x 260 x 250	5,5	63	240 x 120 x 70	10 x 100	
1172-D	365 x 260 x 250	5,8	63	240 x 120 x 70	10 x 100	

EAN Codes

Item	Inner case	Outer case
1172-A	8718007911411	8718007011418
1172-B	8718007911428	8718007011425
1172-C	8718007911435	8718007011432
1172-D	8718007911442	8718007011449