

## 1174 SAFETOUCH® ADVANCED™ PLATINIUM NITRILE POWDER-FREE Class I, Medical Device and Category II, Personal Protective Equipment

### FEATURES & BENEFITS

- Great elasticity
- Good dexterity
- Its smooth appearance gives it correct touch and Grip of objects
- Made from a thinner nitrile material, Slim Nitrile technology provides great protection and flexible “latex-like” comfort
- Watertight and airtight




### Recommendation:

- Glove for professional use
- In the healthcare and industrial environment, this glove is used to make barrier against microorganisms between the wearer of the glove and the examined people or objects / materials handled.
- Powder-free, it reduces cross contamination and limits the risks of irritation.
- For food handling / preparations (HACCP)

### PRODUCT DESCRIPTION

Product name :	<b>SafeTouch® Advanced™ Platinum</b>
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	White
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Textured fingertips
Powdered :	No
Length :	240 mm
Palm thickness :	0,08 mm
Finger thickness :	0,12 mm
Unit weight (M):	4 g
Country of origin :	Thailand
AQL :	< 1,5
Strength at break:	
✓ Before aging	≥6N
✓ After aging	≥6N

Item	Size	
1174-A	T 6	XS
1174-B	T 7	S
1174-C	T 8	M
1174-D	T 9	L
1174-E	T 10	XL



### INDUSTRIES

Health / Medical / Dentist / Tattoo / Food / Non residential catering / Community / Administration / Sanitary / Maintenance / Cleaning / Lab / Pharmaceutical / Research

## CERTIFICATION

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC

Comply with the EN 455-1/2/3/4 Medical gloves for single use

Manufactured under ISO 13485 certified quality system

In accordance with requirements of personal protective equipment Dir.89/686/EEC

Comply with the EN 420 & EN 374-2:2003 : Protection against microorganisms.

In accordance with requirements of Reg.(EC) 1935/2004 and decree of 9th November, 1994. This glove is suitable for food contact\*.

\*Restriction : not put in contact with very fatty foods (see Annex 1).

## STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

Shelf life of the product: 5 years

## LOGISTIC INFORMATION








### Specifications outer case



### Specifications inner case

### Colour code inner and outer

Item	Size mm	Gross weight kg	QTY/pallet	Size mm	QTY	Colour
1174-A	330 x 255 x 228	4,1	60	220 x 125 x 63	10 x 100	
1174-B	330 x 255 x 228	4,7	60	220 x 125 x 63	10 x 100	
1174-C	330 x 255 x 228	5,1	60	220 x 125 x 63	10 x 100	
1174-D	330 x 255 x 228	5,3	60	220 x 125 x 63	10 x 100	
1174-E	330 x 255 x 228	6,2	60	220 x 125 x 63	10 x 100	

### EAN Codes

Item	Inner case	Outer case
1174-A	8 718 007 902 259	8 718 007 002 256
1174-B	8 718 007 902 266	8 718 007 002 263
1174-C	8 718 007 902 273	8 718 007 002 270
1174-D	8 718 007 902 280	8 718 007 002 287
1174-E	8 718 007 902 297	8 718 007 002 294