

1191 MEDICOM® ADVANCED VITALS™ NITRILE POWDER-FREE Medical Device, Class I

FEATURES & BENEFITS

- Latex-like feel
- Agreeable dexterity
- Excellent grip in both wet and dry conditions
- NEW convenient packaging that holds 200 gloves per box reducing storage space, transportation costs and packaging waste



Recommendation:

- In the healthcare and industrial environment, this glove is used to make barrier against microorganisms between the wearer of the glove and the examined people or objects / materials handled.
- Powder-free, it reduces cross contamination and limits the risk of irritation.
- For food handling / preparations (HACCP)

PRODUCT DESCRIPTION

Product name :	Medicom® Advanced Vitals™
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	Blue
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Textured fingertips
Powdered :	No
Length :	240 mm
Palm Thickness :	0,07 mm
Finger Thickness :	0,10 mm
Unit weight (M):	3 g
Country of origin :	Thailand
AQL :	<1,5
Strength at break :	
Before aging	≥ 6N
After aging	≥ 6N

Item	Size	
1191-A	T 6	XS
1191-B	T 7	S
1191-C	T 8	M
1191-D	T 9	L
1191-E	T 10	XL



INDUSTRIES

Food / Non residential catering / Health / Medical / Dental clinic / Community / Administration / Sanitary (except household products)

CERTIFICATION

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC
 Comply with the EN 455-1/2/3/4 Medical gloves for single use
 Manufactured under ISO 13485 certified quality system

In accordance with requirements of Reg.(EC) 1935/2004 and decree of 9th November, 1994. This glove is suitable for food contact*.

*Restriction : not put in contact with very fatty foods (see Annex 1).

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be exposed to moisture and sunlight, must be stored at a temperature between 5°C and 40°C
 Shelf life of the product: 5 years

LOGISTIC INFORMATION








Specifications outer case



Specifications inner case

Colour code inner and outer

Item	Size mm	Gross weight kg	QTY/Pallet	Size mm	QTY	Colour
1191-A	420 x 270 x 260	6,23	49	220 x 120 x 85	10 x 200	
1191-B	420 x 270 x 260	6,54	49	220 x 120 x 85	10 x 200	
1191-C	420 x 270 x 260	7,16	49	220 x 120 x 85	10 x 200	
1191-D	420 x 270 x 260	8,48	49	220 x 120 x 85	10 x 200	
1191-E	420 x 270 x 260	10,35	49	220 x 120 x 85	10 x 200	



EAN Codes

Item	Inner case	Outer case
1191-A	3 700 381 508 795	3 700 381 508 634
1191-B	3 700 381 508 801	3 700 381 508 641
1191-C	3 700 381 508 818	3 700 381 508 658
1191-D	3 700 381 508 825	3 700 381 508 665
1191-E	3 700 381 508 832	3 700 381 508 672