



## HAPPYcover<sup>®</sup>



### SGAIEL POWDER

TECHNICAL DATA SHEET PAGE 1:2

<b>TRADING NAME</b>	<b>HAPPY COVER<sup>®</sup> SGAIEL POWDER</b>
<b>DESCRIPTION</b>	Super absorbent powder cross-linked polyacrylate based.
<b>ITEMS CODES</b>	<p>AM 005 (200 gr)</p> <p>AM 004 (400 gr)</p> <p>AM 006 (1000 gr)</p> <p>AM 010 (5000 gr)</p>
<b>RDM (Italy)</b>	<p>1572422 (AM005 - 200 gr)</p> <p>1572425 (AM004 - 400 gr)</p> <p>1572504 (AM006 - 1000 gr)</p> <p>1638004 (AM010 - 5000 gr)</p>
<b>CND (Italy)</b>	A060399
<b>HS CODE</b>	3906.90.90.00
<b>PACKAGING</b>	1 jar
<b>EXPIRING TIME</b>	5 years
<b>BELONGING CLASS OF MD</b>	Class I medical device - Not sterile Pursuant to 93/42/EEC
<b>CE MARK</b>	Pursuant to 93/42/EEC and 2007/47 EEC of 09/02/2007




## HAPPYcover®



### SGAIEL POWDER

#### TECHNICAL DATA SHEET PAGE 2:2

<b>INTENDED USE</b>	Thickening powder for preventing of liquid spills.
<b>RAW MATERIALS</b>	White cross-linked sodium polyacrylate.
<b>INSTRUCTIONS FOR USE</b>	Pour the powder on spills and wait for the liquids to be thickened. After that remove it and dispose as per current legislation.
<b>CAUTION</b>	Do not use sharp objects to open the package Do not disperse in the environment after use. Do not use if package is open or damaged. Not for human consumption.
<b>CONTRAINDICATIONS</b>	No contraindications by using these medical devices.
<b>STORAGE</b>	Store in a cool, dry place away from heat sources and sunlight..
<b>TRANSPORTATION</b>	Not regulated by EEC and Italia Transport Department as dangerous product. Terrestrial transport (RID/ADR): not dangerous. Maritime transport (IMG): not dangerous. Airplane transport (ICAD/IATA): not dangerous.
<b>APPLICABLE REGULATION</b>	93/42/EEC of 06/14/1993 2007/47 EEC of 09/05/2007
<b>GMDN</b>	Code GMDN: 3906.90.90.00 Description: Powder super absorbent sodium polyacrylate cross-linked based enclosed in water-soluble film
	Angimed S.r.l. Corso Umberto I, 590 ed.66 - 65015 Montesilvano (PE)

CE