

Medicom® Technical Data Sheet

2020 MEDICOM® SAFE+MASK® Cone Mask with Earloop, Blue

PRODUCT DESCRIPTION

Molded contour mask with nosepiece and stretchable headband.



PRODUCT DETAILS

Product name: **Safe+Mask®**
Product type: Single use, non-sterile
Material: Polyester
Nosepiece: Flat Aluminium
Colour: Blue
Quantity per case: 12 x 50 units
Country of origin: China

Item	Colour
2020	Blue



STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun
Storage temperature: ≤ 40°C

N.B.

This mask doesn't meet neither the requirements for medical mask (EN 14683) nor the ones for PPE mask (EN 149), which is why it can not display the CE marking.



Bld de la chanterie | BP 10059 – 49181 Saint Barthélemy d'Anjou Cedex
T : 33 (0) 241 963 434 | F : +33 (0) 241 963 453 | sales@kolmi-hopen.com | www.h-open.com