

2025

SAFE+MASK® PREMIER PROSHIELD™ TYPE IIR – Earloop Medical Device, Class I


FEATURES & BENEFITS

- 3-ply non-woven mask
- Type IIR protective splash-resistant material (according to ISO 22609:2004)
- Polyester visor allows eyes protection against the risk of potentially contaminated liquids projection
- Adjustable nosepiece to customize your own desired fit
- Comfortable, soft, latex-free earloops.



PRODUCT DETAILS

Product name:	Safe+Mask® Premier ProShield™ Type IIR Earloop
Product type:	Single use, non-sterile
Outer/inner layer:	Spunbond polypropylene
Filter:	Meltblown polypropylene
Option:	Visor
Colour:	Teal
Origin:	China

Item	Colour
2025	Teal 



CERTIFICATION

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC.

Manufactured under ISO 13485 certified quality system.

In conformity with the relevant harmonized standard EN 14683 and EN ISO 14971.

Mask performance NF EN 14683:2014		
TEST	Type IIR mask to protect the patient as well as the healthcare provider against splashes of potentially contaminated liquids.	
	Standard	2025
Bacterial Filtration Efficiency (BFE) (%)	≥ 98	≥ 98
Differential pressure (Pa/cm ²)	< 49	< 49
Splash resistance pressure (kPa)	≥ 16	≥ 16

INDUSTRIES

Health / Medical / Surgical / Dental

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

Shelf life of the product: 3 years

LOGISTIC INFORMATION



Specifications outer case

Specifications inner case

Item	Size mm	Gross weight kg	QTY/ pallet	Size mm	QTY
2025	390 x 195 x 235	1,42	105	221 x 187 x 95	4 boxes of 25 pieces



EAN Codes

Item	Outer case	Inner case
2025	8 718 007 003 024	8 718 007 903 027