

Technical Data Sheet

2119BV

TYPE IIR MEDICAL MASK WITH EARLOOPS
Medical Device Class I

PRODUCT DESCRIPTION

Product name: **Mask#53A**
 Product type: Single use, non sterile
 Inner layer: Spunbond polypropylène
 Filter: Meltblown polypropylène
 Outer layer: Spunbond polypropylène
 Links : Earloops
 Dimension : 95 mm (± 5 mm) x 175 mm (± 7 mm)
 Packaging : 10 x dispenser boxes of 50 masks
 Origin : China

Item	Color
2119BV	White <input type="checkbox"/>



Filtration and Comfort Information- NF EN 14683:2019 Type IIR

TEST	REQUIRED LEVEL	LABORATORY	REPORT NUMBER	DATE OF REPORT	AVERAGE	MINIMUM VALUE	MAXIMUM VALUE
BACTERIAL FILTRATION EFFICIENCY: BFE (%)	≥ 98	NELSON	1320119-S01	18/08/2020	99,5	99,3	99,7
DELTA P (Pa/cm ²)	< 60	NELSON	1320119-S01	18/08/2020	47,4	44,9	49,3
SPLASH RESISTANCE (kPa)	≥ 16	NELSON	1320120-S01	06/08/2020	PASS		

PRECAUTIONS

The mask must be used on intact skin only.
 Re-use or extended use may lead to infection or cross-contamination.
 Please refer to your national legislation on the subject of disposal of this product once it has been used.

CERTIFICATION & STANDARDS

Complies with the requirements of Regulation (EU) 2017/745 relating to Medical Devices.
 Complies with applicable harmonized standards EN 14683.
 Manufactured under an ISO 13485 certified system.

STORAGE

Normal conditions of conservation and storage: must not be exposed to humidity and sun, must be stored at a temperature between 5°C and 40°C.
 Product life: 3 years

2119BV

TYPE IIR MEDICAL MASK WITH EARLOOPS

Medical Device Class I

LOGISTICS



Outer case

Inner box

Item	Size mm	Gross weight	QTY/pallet	Size mm	QTY
2119BV	495 x 210 x 215	2 kg	72	103 x 103 x 215	10 boxes of 50 pieces

EAN codes

Item	Outer case	Dispenser box
2119BV	8 718 007 913 903	8 718 007 913 910

www.biomark.nl

