

**2313, 2315, 2316
2317, 2318, 2319**
SAFE+MASK® ECONOMY™ / Masks with earloops

FEATURES & BENEFITS

- 3-ply non-woven mask
- Single pleat design
- Adjustable nosepiece to customize your own desired fit.
- Comfortable, soft, latex-free earloops
- Good price / quality ratio

Recommended use:







- Short, low risk procedures
- It protects the working environment and the products manufactured by limiting the propagation of accidental mouth and nasal sprays (sneezing, coughing ...)



PRODUCT DETAILS

Product name:	Safe+Mask® Economy™
Product type:	Single use, non-sterile
Outer/inner layer:	Spunbond polypropylene
Filter:	Meltblown polypropylene
Link:	Earloops
Size:	95 x 175 mm (± 5 mm)
Unit weight:	2,6 g (± 10%)
Packaging:	Box of 50 units
Origin:	China



Item	2313	2315	2316	2317	2318	2319
Colour	Lavender	Blue	Pink	Yellow	Green	White
						

INDUSTRIES

Petrochemical / Food / Away from home food service / Non residential catering / Community / Administration / Hygiene / Maintenance / Cleaning / Lab / Pharmaceutical / Research

2313, 2315, 2316
2317, 2318, 2319

SAFE+MASK® ECONOMY™ / Masks with earloops

CERTIFICATION & STANDARDS

Manufactured under EN ISO 9001 certified quality system

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

Shelf life of the product: 3 years

LOGISTICAL INFORMATION



Outer case specifications



Inner case specifications

Item	Size mm	Gross weight kg	QTY/pallet	Size mm	QTY
2313 2315 2316 2317 2318 2319	510 x 210 x 165	1,9	80 cartons (10 layers of 8 cartons)	190 x 100 x 80	10 boxes of 50 pieces

EAN Codes

Item	Outer case	Inner case
2313	8 718 007 910 711	8 718 007 910 704
2315	8 718 007 012 781	8 718 007 912 784
2316	8 718 007 910 735	8 718 007 910 728
2317	8 718 007 910 759	8 718 007 910 742
2318	8 718 007 912 798	8 718 007 912 791
2319	8 718 007 910 773	8 718 007 910 766