



## 2950

### FLOWER MASKS Type II with earloops Medical Device, Class I

#### FEATURES & BENEFITS

- Protects the environment from droplets emitted by the wearer.
- High filtration with filter media (polypropylene).
- 2 comfortable large earloops.
- Very soft nonwoven.
- Anatomic folding limiting inward leakage.
- Integrated nose piece.



#### RECOMMENDED USE

- Intended to be used by healthcare professionals.
- Medical face masks may also be intended to be worn by patients and others to reduce the risk of spreading infection, particularly in epidemic or pandemic situations.

#### PRODUCT DETAILS

Product name: **Flower Masks - Type II with earloops**  
 Product type: Single use, non-sterile  
 Straps: Elastic Earloops  
 Size: 90 x 175 mm (± 5 mm)  
 Unit weight: 2,8 g (± 10%)  
 Packaging: Box of 60 units  
 Country of origin: France



<b>Item</b>	2950-DB	2950-LM	2950-LV	2950-PK	2950-TU
	Dark blue	Lime	Lavender	Pink	Turquoise

**Colour**



#### Product Performance- NF EN 14683:2019 + AC:2019 Type II

Measurements done by NELSON Laboratory (report available on request)

TEST	Required level	Laboratory	Report number	Date of Report	Average	Minimum value	Maximum value
Bacterial Filtration Efficiency: BFE (%)	≥98	NELSON	1250735-S01	07/01/2020	98,8	98,6	99,2
DELTA P (Pa/cm²)	<40	NELSON	1250735-S01	07/01/2020	23,4	22,6	24,3

# 2950

## FLOWER MASKS Type II with earloops Medical Device, Class I

### CERTIFICATION & STANDARDS

Meet the requirements of Regulation (EU) 2017/745 on medical devices.  
Comply with the applicable standards EN 14683:2019 + AC: 2019.  
Manufacturing site: ISO 9001 and ISO 13485 certification.

### PRECAUTIONS

The mask must be used on intact skin only.  
The reuse or extended use may lead to infection or cross-contamination.  
Please refer to your national legislation on the subject of disposal of this product once it has been used.

### STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored at a temperature between 5°C and 40°C.  
Shelf life of the product: 5 years.

Marking on Packaging:



\* Do not use if package is damaged

### LOGISTICAL INFORMATION

Outer case specifications

Inner case specifications

Item	Size mm	Gross weight kg	QTY/pallet	Size mm	QTY
2950-DB 2950-LM 2950-LV 2950-PK 2950-TU	400 x 400 x 220	3,6	48 cartons (8 layers of 6 cartons)	211 x 97 x 97	16 boxes of 60 pieces

### EAN Codes

Item	Outer case	Inner case
2950-DB	8 718 007 910 223	8 718 007 910 230
2950-LM	8 718 007 910 285	8 718 007 910 292
2950-LV	8 718 007 910 247	8 718 007 910 254
2950-PK	8 718 007 910 261	8 718 007 910 278
2950-TU	8 718 007 910 209	8 718 007 910 216

