

EU Type-Examination Certificate

Certificate number: 2777/11723-01/E05-01

This EU Type-Examination Certificate covers the following product group(s) supported by testing to the relevant standards/technical specifications and examination of the technical file documentation:

Following the EU Type-Examination this product group has been shown to satisfy the applicable essential health and safety requirements of Annex II of the PPE Regulation (EU) 2016/425 as a Category III product.

Product reference:

Nitril LIOX

Description:

Five fingered, chlorinated, antimicrobial powder free nitrile gloves, with textured fingers and beaded cuff. Thickness \geq 2.2mil

Sizes:

- 5 – 5.5 (XS)
- 6 – 6.5 (S)
- 7 – 7.5 (M)
- 8 – 8.5 (L)
- 9 – 9.5 (XL)

Colour: Blue

Classification:

EN ISO 374-1:2016 Type C

- Methanol (A)
- 40% Sodium Hydroxide (K)
- 25% Ammonium Hydroxide (O)
- 30% Hydrogen Peroxide (P)
- 37% Formaldehyde (T)

Permeation Level

- 0
- 6
- 0
- 1
- 4

EN 374-4:2013 Degradation (%)

- 97.8
- 2.8
- 19.8
- 42.3
- 25.5

EN ISO 374-5:2016

- Protection against bacteria and fungi
- Protection against viruses

Level

- Pass
- Pass

Standards/Technical specifications applied:

EN 420: 2003+A1: 2009; EN ISO 374-1:2016+A1:2018; EN ISO 374-5:2016; EN 421: 2010 (excluding clause 4.3)

Technical reports/Approval documents:

CHM0262245/1738/EN/B, CHM0264616/1747/CL/A, CHM0262245/1738/EN/A, CHM0264616/1747/SPT, CHM0264616/1747/CL/B, CHM0274345/1833/LC/B

Signed on behalf of SATRA:



Quincey Brown



Geoff Graham

Date of issue: 26/02/2019

Expiry date: 20/12/2023

TERMS AND CONDITIONS

The following conditions apply in addition to SATRA's standard terms and conditions of business and those given in the current certification agreement.

The certificate holder is licensed to mark the products detailed within this certificate in accordance with Annex V (Module B) of the Regulation (EU) 2016/425 of the European Parliament and of the council of 9th March 2016 on personal protective equipment once you have drawn up an EU declaration of product conformity. Please note:

1. Where the product is classified as category III then CE Marking of production is reliant on current compliance with Regulation 2016/425 module C2 or Module D. (Except that specifically produced to fit an individual user).
2. Full details of the certification and product are contained within the manufacturer's technical documentation.
3. Where a translation of this certificate exists, the English language version shall be considered as the authoritative text.
4. Certification is limited to production undertaken at the sites listed in the manufacturers technical documentation.
5. Ongoing manufactured product shall be consistent with the product(s) certified and listed on this certificate.
6. The Manufacturer shall inform SATRA of any changes to the certified product or technical documentation.
7. This certificate shall be kept together with the relevant technical documentation in a safe place by the client named on this certificate. Production of this certificate and other documentation may be required by a representative of the EC member state government.
8. This certificate relates only to the condition of the testable items at the time of the certification procedure and is subject to the expiry date shown.
9. SATRA Technology reserves the right to withdraw this certificate if it is found that a condition of manufacture, design, materials or packaging have been changed and therefore no longer comply with the requirements of Regulation 2016/425.