




TECHNICAL DATA SHEET
COUCH ROLLS
Paper+hdpe



NAME OF THE PRODUCT: DISPO SHEET

DESCRIPTION OF THE PRODUCT: Drapes in roll

DIMENSION: H50X50mt, H50x100mt, H60x50mt,
H60x100mt precut 38cm (tolerance +/-5%)

CLASSIFICATION: Medical Device Class I
Classification according to the Annex IX paragraph III rule 1
of the European Directive 93/42/EEC (and successive amendments
included Directive 2007/47/EEC)

TECHNICAL FEATURES: Pe-coated, precut examination drapes in roll, made of one ply of very absorbent tissue and one ply of waterproof polyethylene. Coupling and embossing of the two materials assure a high tearing strength and a perfect waterproofness Latex free
Paper (100% cellulose) 26 gr/m² (+/- 5%)
Polyethylene 7/8 μ/m² (+/- 8%)

MANUFACTURER: Dispotech srl Via M. del Grosso 8-12
23022 Chiavenna (SO)
Tel. 0039.0343.36711
Fax 0039.0343.3656
e-mail: info@dispotech.it

COLOR: White

LOT: Printed on the label applied on the box

STORAGE: Keep at ambient temperature, in a dry place, away from heat sources, flames or sparks.

PACKAGING: each roll wrapped in PE bag, 6 rolls in each box

DATE: 23/11/2018
Rev.03
MADE BY RGQ
APPROVED BY DG

Dispotech srl

Sede operativa: 23020 GORDONA (SO) Via al piano, 29 - Tel. +39 0343 36711 - Fax. +39 0343 36567 - www.dispotech.com
Sede legale: 23022 CHIAVENNA (SO) Via M. Del Grosso, 8/12 - e-mail: info@dispotech.it - PEC: dispotech@legalmail.it
Cap Soc. € 1.500.000,00 i.v. - P. IVA 00672170149 - R.E.A. 47213 C.C.I.A.A. di SO - Uff. Reg. Imp. SO 00672170149